

# DataPaq<sup>™</sup> EXPRESS-NMR Rack Scanner

Fast and compact – scans all 2D-barcoded NMR tubes

### FEATURES

- Industry-leading fast scan and decode time for a full rack of under 3 seconds
- Fast camera-based scanner
- Can be integrated with large robotic systems
- Available with either USB or RJ45 connector
- Capable of scanning all 2D-barcoded NMR tube caps
- Patented technology allows imaging under difficult lighting conditions
- DP5 or DataPaq<sup>™</sup> software with free lifetime upgrades
- Easy to set-up and easy to use

#### Industry-leading fast scan and decode time

Using the smallest footprint 2D barcode rack scanner on the market enabled the unique design of the Ziath Express NMR tube reader. The DataPaq<sup>™</sup> Express 2D scanner offers easy, out-ofthe-box set up and is delivered pre-calibrated and ready to scan. This instrument can scan all 2D-barcoded NMR tubes including Bruker<sup>™</sup> and Norell<sup>™</sup> and is supplied with a handy slide-out tray designed to hold Bruker tube racks in the correct orientation and in the correct position for successful barcode scanning.

# SPECIFICATIONS

- Imager & Housing = 170mm x 210mm x 230mm (WxDxH)
- Controller = 95mm x 195mm x 55mm (WxDxH)
- Total Weight = 2.4kg
- Power adaptor input = 100-240VAC; output = 5VDC, 4A
- Compatible with Bruker<sup>™</sup>, Wilmed<sup>™</sup> and Norell<sup>™</sup> 2D barcoded NMR tubes
- Export data to Excel, text, JSON, XML and image files
- Operating system on host PC = Windows<sup>®</sup> 7, 8, 10

- Data output connector = USB3 or RJ45
- Warranty period = 2 years
- Software includes remote control capability for integration applications



NMR is now an established and vital part of determining the structure of novel compounds in both Pharma and Agri-bio. The task of keeping track of large numbers of similar-looking NMR tubes is now easier as manufacturers now produce NMR tubes with 2D barcodes on their caps, allowing easy identification and tracking of individual tubes. However, due to the barcodes being placed on the top of the tubes, these can be difficult to read with a conventional scanner designed to read from underneath. In addition, there are at least two manufacturers using different formats and colours for their caps and barcodes.

To overcome this Ziath, world leader in 2D bar-code reading technology, have a solution that puts the Ziath Express scanner on top of the tubes without applying any weight, which could damage them. It also provides a dedicated slide for Bruker tube carriers which can only be inserted in one orientation – conserving the A1 position each time. Ziath software has templates for most types of NMR tube and can switch easily between. The Express-NMR is compatible with Datapaq 3.17 or the newer web-based DP5 control software. The Express controller is conveniently separated so can be hidden away under bench without taking up valuable space.

## **INCLUDED WITH THE DataPaq™ EXPRESS-NMR**

- DataPaq<sup>™</sup> Express instrument with separate controller module
- Mounting housing for inverted Express optical module
- Precision engineered slide and drawer with dedicated rack holder
- Power adaptor and all necessary connecting cables

Full data processing options including export of the tube barcodes in tables or as Text/.CSV, Image, XML, XLS and in JSON formats is supported in the easy-to-use Datapaq or DP5 software supplied with the instrument.

### **Ordering information**

Code	Description	Quantity
ZTS-EXP-NMR	DataPaq™ Express-NMR Rack Scanner	1



**Ziath Ltd**, Unit 2A Solopark Trading Estate, Pampisford, Cambridge, CB22 3HB, UK Telephone: +44 (0) 1223 855021 Email: info@ziath.com www.ziath.com