



# **DATAPAQ™ SINGLE TUBE**

Small and cost effective 2D  
barcode scanner



## **FEATURES**

- Fast scanning of 2D tubes
- Compatible with all 2D Datamatrix tubes
- 'Keyboard Wedge' setting allows data entry directly into your application
- Strong mineral glass window
- DataPaq™ software with free lifetime upgrades
- Cryoprotection option

## **Easy to set-up and easy to use**

- Fast scanning of 2D tubes
- Compatible with all 2D Datamatrix tubes
- 'Keyboard Wedge' setting allows data entry directly into your application
- Strong mineral glass window
- DataPaq™ software with free lifetime upgrades
- Cryoprotection option
- The DataPaq™ Single Tube 2D scanner offers easy, out-of-the-box set up and is delivered pre-calibrated and ready to scan
- This instrument can scan all commercially available 2D datamatrix tubes
- The Single Tube scanner can be used for applications that require data input without direct interaction with your computer

# SPECIFICATIONS

- Dimensions = 60mm x 90mm x 60mm (WxDxH)
- Total Weight = 0.4kg
- Powered by USB - connect directly into computer
- Operating system on host PC = Windows® 7, 8, 10
- Warranty period = 2 years
- Free DataPaq™ software upgrades for life
- Software includes remote control capability for integration applications

## INCLUDED WITH THE SINGLE TUBE

- DataPaq™ single tube instrument
- DataPaq™ single software
- User guide



### Option

#### CRYOPROTECTION

Prevents the scanner window from misting up when a cold tube is scanned. Our passive technology does not use heat or blown air.

## Ordering information

Code	Description	Quantity
ZTS-S2	DataPaq Single Tube Scanner	1
ZTS-S2-80	DataPaq Single Tube Scanner with Cryoprotection	1



**Ziath Ltd**, Unit 2A Solopark Trading Estate, Pampisford, Cambridge, CB22 3HB, UK  
Telephone: +44 (0) 1223 855021 Email: [info@ziath.com](mailto:info@ziath.com) [www.ziath.com](http://www.ziath.com)