



PRESS RELEASE

Compact Rack Scanner for Easy Liquid Handling Robot Integration

Compatible with all 2D barcoded tube racks in SBS format, the **Ziath DataPaq™ Express** is the smallest footprint rack scanner currently available.

With a footprint of only slightly more than one plate / rack position, the DataPaq™ Express offers a **space-saving way to integrate a fast full-rack scanner on to a liquid handling robot**. Designed with a separate power and processing box, which can be positioned under your liquid handling robot and under your deck, using a DataPaq™ Express scanner frees up vital deck space.

The uniquely low form factor of the DataPaq™ Express allows liquid handling robots easy gripper access to simply pick up and dispense from racks on top of the scanner. Offering rapid image scanning and decoding in just two seconds, the camera based DataPaq™ Express will also help improve your robotic workflow.

Baseplates and drivers are available to enable easy integration of the DataPaq™ Express with most commercial liquid handling robots, including those from Tecan, Hamilton, Perkin Elmer, Beckman, Andrew Alliance and Hi-Res Bio.

Ziath Ltd
Unit 2a
Solopark Trading Estate
Station Road
Pampisford
Cambridge
CB22 3HB
United Kingdom
Tel: +44 (0)1223 855021
Email: queries@ziath.com
Web: www.ziath.com



The DataPaq™ Express scanner is easy to set-up and can be ready to use only a few minutes after taking it out of the box. The fast full-rack scanner arrives fully calibrated and is ready to read all makes of SBS format 2D barcoded racks and tubes. Powerful DataPaq™ software makes it simple to export data to Excel, XML or text/CSV, and scanned images can also be saved. The DataPaq™ software can also connect with Oracle, SQL Server, MySQL, Postgres and other databases.

For further information please visit <http://www.ziath.com/index.php/products/barcode-scanning/datapaq-express-rack-reader> or contact Ziath on +1- 858-880-6920 / +44-1223-855021 / info@ziath.com.

Founded in 2005, Ziath specialises in development of innovative instrumentation control and information management products using 2D DataMatrix bar-coded tubes to simplify automation processes in life science organisations, from academia, to the biotech and pharma industries. For further information please visit www.ziath.com.


DECEMBER 2020

ziathpr25

Ziath Ltd
Unit 2a
Solopark Trading Estate
Station Road
Pampisford
Cambridge
CB22 3HB
United Kingdom
Tel: +44 (0)1223 855021
Email: queries@ziath.com
Web: www.ziath.com



Illustrative images (available on request)



For further information:

MEDIA: Dr Bill Bradbury
TECHNICAL: Steve Knight

Primetek Solutions
Ziath Ltd

tel. +44-208-546-0869 /
tel. +44-1223-855021 /

email info@primetek-solutions.com
email info@ziath.com

Ziath Ltd. Registered Office: Unit 2a Solopark Trading Estate, Station Road, Pampisford, Cambridge, CB22 3HB, UK

Registered in the UK (no. 05598930) VAT - GB890276990